

Medi-Table

Prep / Surgery Table



Operating / Prep Table Hydraulic

The ergonomic benefits of variable height prep/surgery tables are perfect for lifting large animals to a comfortable working operating height.

The Medi-Table VT200 (Hydraulic) Series and VT201 (Electric) allows elevation by hydraulic or electric actuator to a convenient working height and includes a tilt mechanism.

- Reliable lifting system (3-Year Warranty)
- Tilt Mechanism
- Max. Working Height 110cm
- Choice of wheels or level foot system with retractable castors
- Tough waterproof melamine top (warmer than SS)
- Choice of working lengths (120, 150, 180cm)
- Width 60cm
- Lifting capacity 260Kg (Electric) 180Kg Hydraulic
- High-density rubber mat

Options

- Battery Power System (charge overnight and enables 20-50 lifting cycles per day) – NO MORE unsightly and nuisance cables during consultation.
- IV Pole Attachment
- Tie-Down Options
- Slide out panel for monitoring equipment