



Choose the right stethoscope to help you perform at your best.

3M™ Littmann® Stethoscopes feature adjustable headsets with soft-sealing eartips. Select Littmann mechanical stethoscopes are designed with tunable technology to give you the ability to hear different frequency sounds by simply adjusting the pressure applied to the chestpiece. Get the acoustic performance, high quality and comfort you need with a Littmann stethoscope.



	Lightweight II S.E. Stethoscope	Classic II Paediatric and Infant Stethoscopes	Master Classic II™ Stethoscope	Classic III™ Monitoring Stethoscope	Cardiology IV™ Diagnostic Stethoscope	Master Cardiology™ Stethoscope	CORE Digital Stethoscope
	Ideal for taking blood pressure readings and making limited patient assessments.	Hard-working, dependable, clinical tools for physical assessment and monitoring of your patients, one long shift after another.			Provides outstanding acoustic sensitivity to help you perform detailed diagnostic auscultation. Useful in challenging and critical care environments.		High-performance auscultation. Hear more with amplification, noise cancellation and mobile connectivity.
Approximate Length	71cm	71cm	69cm	69cm	69cm	69cm	69cm
Approximate Weight	118g	105g, 95g	160g	150g	177g	185g	232g
Warranty	2 years	3 years	3 years	5 years	7 years	7 years	2 years
Diaphragm	Tunable	Non-Tunable (Floating)	Tunable	Tunable (2)	Tunable (2)	Tunable	Tunable (2)
Approximate Chestpiece Diameter	Double-sided 5.2 x 4.5cm 3.5 x 3.2cm	Double-sided Paediatric 3.4cm 2.7cm Infant 2.8cm 1.7cm	Single-sided 4.3cm	Double-sided Adult side 4.3cm Paediatric side 3.4cm	Double-sided Adult side 4.3cm Paediatric side 3.4cm	Single-sided 4.8cm	Double-sided Adult side 4.3cm Paediatric side 3.4cm
Features	Teardrop-shaped chestpiece Fits under blood pressure cuff	Specially designed for the smallest patients	Designed for easier grip and maneuvering	Convertible open bell	Convertible open bell Dual-lumen tubing	Dual-lumen tubing Best acoustic response in Littmann mechanical line	Up to 40x amplification* Switch from analog listening to digital Connects to Eko software on a smart device** to capture, store and share data



The power to learn,
in the palm of your hand.

3M™ Littmann® Learning Institute App –
Free with purchase of a 3M™ Littmann® Stethoscope!

Need more help? Visit Littmann.co.uk

*At peak frequency, vs. analogue mode
**Not included



3M™ Littmann® Stethoscope Ordering Information and Tube Colour Guide

Components
(Chestpiece/Eartube/Stem)

Lightweight II SE Stethoscope	Classic II Infant Stethoscope	Classic II Paediatric Stethoscope	Master Classic II™ Stethoscope	Classic III™ Monitoring Stethoscope	Cardiology IV™ Diagnostic Stethoscope	Master Cardiology™ Stethoscope	CORE Digital Stethoscope
-------------------------------	-------------------------------	-----------------------------------	--------------------------------	-------------------------------------	---------------------------------------	--------------------------------	--------------------------

Tube Color

All Standard Finish	2450 Black	2114 Black	2113 Black	2144L Black	5620 Black	5627 Burgundy	6152 Black	6156 Plum	2160 Black	
	2451 Burgundy	2114R Red	2113R Red	2146 Burgundy	5621 Gray	5633 Pearl Pink	6154 Navy	6158 Raspberry	2163 Burgundy	
	2452 Caribbean	2124 Caribbean	2119 Caribbean	2147 Navy	5622 Navy	5831 Plum	6155 Hunter Green	6184 Burgundy	2164 Navy	
	2454 Ceil Blue		2122 Raspberry		5623 Caribbean	5832 Lavender			2167 Plum	
	2456 Pearl Pink				5630 Ceil Blue	5835 Turquoise				
				5648 Raspberry	5839 Lemon-Lime					
All Black			2141 Black	5803 Black	5868 Burgundy	6163 Black	6168 Navy	2161 Black	8490 Black	
All Brass Finish								2175 Black		
All Copper Finish					5809 Chocolate					
All Smoke					5811 Black			2176 Black		
Rainbow/Brass/Brass	2157 Raspberry	2153 Caribbean			5806 Raspberry	5807 Caribbean				
Brass/Smoke/Smoke										
Rainbow/Black/Black					5870 Black		6165 Black			
Mirror/Smoke/Smoke						5863 Navy				
Mirror/Standard/Mirror					5959 Ceil Blue					
Champagne/Smoke/Smoke					5861 Black	5864 Burgundy	6177 Black	6170 Burgundy		
Champagne/Black/Champagne							6176 Burgundy	6190 Caribbean		
High-Polish Rainbow/Black/Black							6242 Navy			
High-Polish Rainbow/Smoke/Smoke							6240 Black	6241 Raspberry		
COLOURED STEM SERIES: (square inset represents stem colour)										
Black/Black/Colour							6200 Black	6203 Black		
High-Polish Smoke/Black/Colour							6201 Black			
High-Polish Smoke/Smoke/Colour							6202 Navy	6204 Black		
Rainbow/Black/Colour							6232 Black	6234 Caribbean		
Champagne/Champagne/Colour							6238 Gray			
Smoke/Smoke/Colour					5872 Turquoise	5875 Lime				
Mirror/Smoke/Colour					5873 Gray					
Mirror/Standard/Colour					5960 Plum	5962 Pearl Pink				

Components:
Chestpiece/Eartube/Stem

Standard Finish



Black



Brass Finish




Copper Finish



Smoke



Rainbow




Mirror



Champagne



High-Polish Smoke



High-Polish Rainbow



Colours shown are approximate.

www.littmann.co.uk | Also available: Veterinary Stethoscope (1392), Teaching Stethoscope (2138)

© 3M 2021. All rights reserved. 3M, Littmann, the L Littmann logo, Classic III, Cardiology IV, Master Classic II and Master Cardiology are trademarks of 3M. Used under license in Canada. Eko trademarks, including CORE, are marks and/or registered marks of Eko Devices Inc. All other trademarks are property of their respective owners. OMG154398.