



## Product Certificate

Jørgen Kruuse A/S hereby declare the following details regarding the mentioned product.

**Catalogue Number**

271715

**Product name**

BUSTER ICU Heat Mat 12 V

**Product type**

Soft and flexible heating pad with heating settings from 35 to 45°C.  
The mattress is waterproof and can be scrubbed with disinfectants.  
Can be set from 1 to 9 hours.

For veterinary use only.

**Materials**

PVC and polyurethane.

**Size**

48 x 80 cm

**Storage condition**

Store in dry areas (5 - 40°C). Avoid direct sunlight.

**Shelf-life**

N/A

**Country of origin**

China

**Last revised**

13-feb-2023

**Version**

1.0

**Signed**

Michella Lund

**Function**

Quality & Regulatory  
Specialist